



Apparatus Mercury Triple Point

Purpose Designed for Isotech Mercury Cell

- Outstanding Convenience and Safety
- Liquid Free

The Isotech Model ITL -M-17725 Cryostat is a selfcontained, mechanically-refrigerated, system with a main well to house one cell and two auxiliary wells for prechilling of thermometers. The cryostat temperature is steplessly adjustable from -36°C to -42°C with scaledivisions of 0.1°C (interpolation possible) and 24 hour stability better than ± 0.05 °C measured by an SPRT in the well of a Mercury Cell inside the Cryostat.

The cryostat has several unique features providing outstanding convenience and safety. The refrigeration system has sufficient capacity to bring a cell to operating temperature in about one hour. At operating temperature, the cooling rate is about 1 Kelvin/minute and the heating rate is about 2 Kelvin/minute.

This permits rapid changes to be imposed on the temperature of the cell environment to a void excessive demands on the (low) heat-of-fusion energy of the mercury within the cell.

In addition, all temperature control is accomplished through control of refrigerant flow, providing inherently fail-safe operation. Indicators provided for the cryostat are "POWER ON" and "COOLING".

The cryostat provides convenient conditions for operating mercury fixed point cells both in heating and in cooling mode.



| Model | ITL-M-17725 Apparatus |
|-------------------|---|
| Temperature Range | -36°C to -42°C |
| Uncertainty | 0.22mK (with cell) |
| Ambient Limits | 18°C to 28°C |
| Plateau Duration | 8-12 hour Plateau |
| Power | 750W typical. 208-240 VAC, 50/60Hz |
| Dimensions | Height 960mm Width 600mm Depth 560mm Weight 96kg |

How to order ITL-M-17725 Mercury Triple Point Apparatus